


# Resolution studies according to flow rate and rotation speed with piperine purification on Analytical FCPC® with 50 ml column

Kromaton application laboratory, Angers, France

## Introduction


Piperine ( $C_{17}H_{19}NO_3$ ) is the alkaloid responsible for the taste and smell of black pepper. It has been also used for the traditional medicine and like insecticide.




Trials were done on Analytical FCPC® equipped with 50ml column to follow resolution of the separation at different flow rate and rotation speed with few 100 mg injection of crude *Piper nigrum* extract. Detection is done with UV/Vis detector at 254 and 300 nm.



## Results



UV 254 nm FCPC® chromatograms of 100 mg *Piper nigrum* extract injection  
Parameters : 4 and 8 ml/mn, 2000 rpm



UV 300 nm FCPC® chromatograms of 100 mg *Piper nigrum* extract injection  
Parameters : 1500 and 2000 rpm, 6 ml/mn

| Parameters        |                   |
|-------------------|-------------------|
| Flow-rate         | 4 and 8 ml/mn     |
| Rotation speed    | 1500 and 2000 rpm |
| Quantity injected | 100 mg            |
| Solvent system    | ARIZONA           |

| Results             |              |
|---------------------|--------------|
| Separation time     | 15 and 30 mn |
| Solvent consumption | 170 mL       |

## Conclusions

Chromatogram shows same resolution between 4 and 8 ml/mn trials and better efficiency at 2000 rpm than 1500 rpm. In this case, separation with ARIZONA solvent system is achieved very fast at 8 ml/mn and 2000 rpm with injection up to 1 gr.